



Product data sheet REF 45.30.20

11.12.2025

General data			
article number	45.30.20	inductive charging device 110-240V f stand hdl w LED	
Unit	piece	Customs code (HS code)	85044085
GMDN	60909	GTIN	04046826534622
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	
Latex free	yes	Preparation instructions	MPK14
 			
Medicon media	453.11.01 (Fiber Optic Laryngoscopes Brochure)		
Complete picture		Detail images	
size	110-240V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		