



## Product data sheet REF 45.30.23

11.12.2025


General data			
article number	45.30.23	Macintosh Flexitip Blades	
Unit	piece	Customs code (HS code)	90189084
GMDN	46827	GTIN	04046826535285
Risk category	1 (EU MDR)	UMDNS-Code	15-076
Classification	10 (EU MDR)	Base UDI-DI	4046826;N004;05;XR
Net weight (kg)	0,000	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	



Medicon media 453.11.01 (Fiber Optic Laryngoscopes Brochure)

## Complete picture

## Detail images



McINTOSH - FlexiTip	
45.30.21	a = 90 mm
45.30.22	a = 110 mm
45.30.23	a = 130 mm

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	1.4301 (X5CrNi18-10)