



## Product data sheet REF 45.30.24

11.12.2025

General data			
article number	45.30.24	Autoclavable LED Larygoskope Handle Medium 2.5V	
Unit	piece	Customs code (HS code)	90189084
GMDN	46830	GTIN	04046826535155
Risk category	1 (EU MDR)	UMDNS-Code	15-076
Classification	10 (EU MDR)	Base UDI-DI	4046826;N001;06;WT
Net weight (kg)	0,000	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	453.11.01 (Fiber Optic Laryngoscopes Brochure) 453.32.00 (FlexiTrac - Abdominal retractor system Brochure)		
Complete picture		Detail images	
size	2.5V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		