



Product data sheet REF 45.30.25

11.12.2025

General data			
article number	45.30.25	Autoclavable LED Larygoskope Handle small 2.5V	
Unit	piece	Customs code (HS code)	90189084
GMDN	46830	GTIN	04046826535186
Risk category	1 (EU MDR)	UMDNS-Code	15-076
Classification	10 (EU MDR)	Base UDI-DI	4046826;N001;06;WT
Net weight (kg)	0,000	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	453.11.01 (Fiber Optic Laryngoscopes Brochure)		
Complete picture		Detail images	
			
size	2.5V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		