



Product data sheet REF 45.30.28

25.08.2026

General data			
article number	45.30.28	Endcap Medium Laryngoscope Handles 2.5V Satin	
Unit	piece	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826534905
Risk category	N.A. (EU MDR)	UMDNS-Code	
Classification	NA (EU MDR)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	
Latex free	yes	Preparation instructions	
			
Medicon media			
Complete picture		Detail images	
size	2.5V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		