



Product data sheet REF 45.30.30

25.08.2026

General data			
article number	45.30.30	Endcap for 3.5V Small Laryngoscope Handles, Satin, Brass C/P. 3.5V Satin	
Unit	piece	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826535254
Risk category	N.A. (EU MDR)	UMDNS-Code	
Classification	NA (EU MDR)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	
Latex free	yes	Preparation instructions	



Medicon media

Complete picture

Detail images

size	3.5V	Coating	
------	------	---------	--

Storage conditions

No special requirements

Materials

Material number (DIN designation)