





## Product data sheet REF 45.48.65

11.12.2025

General data			
article number	45.48.65	replacement lamp 2.7V small p/6	
Unit	box/6	Customs code (HS code)	85392198
GMDN	60909	GTIN	04046826048594
Risk category	KEIN MP (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,012	Operating instructions	
Latex free	yes	Preparation instructions	KEIN MP
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT )		
Complete picture		Detail images	
			
size	2.7V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		