





Product data sheet REF 45.48.66

11.12.2025

General data			
article number	45.48.66	replacement lamp 2.7V large p/6	
Unit	box/6	Customs code (HS code)	85392198
GMDN	60909	GTIN	04046826048600
Risk category	KEIN MP (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,014	Operating instructions	
Latex free	yes	Preparation instructions	KEIN MP
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT)		
Complete picture		Detail images	
			
size	2.7V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		