



Product data sheet REF 45.48.70

11.12.2025

General data			
article number	45.48.70	krypton lamp 2.5V only	
Unit	piece	Customs code (HS code)	85392198
GMDN	60909	GTIN	04046826048617
Risk category	KEIN MP (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,002	Operating instructions	
Latex free	yes	Preparation instructions	KEIN MP
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT)		
Complete picture		Detail images	
			
size	2.5V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		