



## Product data sheet REF 45.52.02

11.12.2025

| General data                                                                                                                                                                  |                                                                                |                           |                    |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|---------------------------|--------------------|
| article number                                                                                                                                                                | 45.52.02                                                                       | Davis Meyer mouth gag cpl |                    |
| Unit                                                                                                                                                                          | piece                                                                          | Customs code (HS code)    | 90189084           |
| GMDN                                                                                                                                                                          | 35085                                                                          | GTIN                      | 04046826049041     |
| Risk category                                                                                                                                                                 | 1 (EU MDR)                                                                     | UMDNS-Code                | 11-822             |
| Classification                                                                                                                                                                | 5 (EU MDR)                                                                     | Base UDI-DI               | 4046826;N004;05;XR |
| Net weight (kg)                                                                                                                                                               | 0,670                                                                          | Operating instructions    | GATDN00X           |
| Latex free                                                                                                                                                                    | yes                                                                            | Preparation instructions  |                    |
| <div></div> |                                                                                |                           |                    |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT ) |                           |                    |
| Complete picture                                                                                                                                                              |                                                                                | Detail images             |                    |
|                                                                                                                                                                               |                                                                                |                           |                    |
| Storage conditions                                                                                                                                                            |                                                                                |                           |                    |
| No special requirements                                                                                                                                                       |                                                                                |                           |                    |
| Materials                                                                                                                                                                     |                                                                                |                           |                    |
|                                                                                                                                                                               | Material number (DIN designation)                                              |                           |                    |
| Main part                                                                                                                                                                     | 1.4021 (X20Cr13)                                                               |                           |                    |
| Handle                                                                                                                                                                        | 1.4301 (X5CrNi18-10)                                                           |                           |                    |
| Others                                                                                                                                                                        | 1.4310 (X10CrNi18-8)                                                           |                           |                    |