






Product data sheet REF 46.37.47

11.12.2025

| General data | | | |
|---|--|----------------------------|--------------------|
| article number | 46.37.47 | needle case 65x42x8mm perf | |
| Unit | piece | Customs code (HS code) | 73269098 |
| GMDN | 12143 | GTIN | 04046826049836 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 16-767 |
| Classification | 1 (EU MDR) | Base UDI-DI | 4046826;N003;03;X8 |
| Net weight (kg) | 0,045 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 65x42x8mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |