






Product data sheet REF 47.27.34

11.12.2025

| General data | | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|--------------------------|--------------------|
| article number | 47.27.34 | metal tray 340x210x60mm | |
| Unit | piece | Customs code (HS code) | 73269098 |
| GMDN | 42893 | GTIN | 04046826050771 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 15-561 |
| Classification | 1 (EU MDR) | Base UDI-DI | 4046826;N003;03;X8 |
| Net weight (kg) | 0,711 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.40 (Katalog Nr 5 Dental) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 340x210x60mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |