




Product data sheet REF 49.40.39

25.08.2026

General data			
article number	49.40.39	Medicontainer 300x135x47mm gold	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826051372
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	0,640	Operating instructions	IFU-0003
Latex free	yes	Preparation instructions	
			
Medicon media			
Complete picture		Detail images	
			
size	300x135x47mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.0255_3.3315 (Al 99,5% + AlMg1)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4301 (X5CrNi18-10)		
Others	SILIKON (Poly(organo)siloxane)		