








Product data sheet REF 49.40.50

11.12.2025

General data			
article number	49.40.50	Medicontainer 300x135x62mm silver	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826051389
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	0,660	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.40 (Katalog Nr 5 Dental) 451.02.61 (brochure Medicontainer System)		
Complete picture		Detail images	
  			
size	300x135x62mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.0255_3.3315 (Al 99,5% + AlMg1)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4301 (X5CrNi18-10)		
Others	SILIKON (Poly(organo)siloxane)		