

Product data sheet REF 49.40.61


25.08.2026

General data			
article number	49.40.61	Medicontainer 300x135x71mm blue	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826051440
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	0,660	Operating instructions	IFU-0003
Latex free	yes	Preparation instructions	



Medicon media

Complete picture	Detail images
------------------	---------------



Seq.	Seq.	Seq.	Seq.	Seq.	4 years	5 years	6 years
400.40.30	400.40.30	400.40.37	400.40.38	400.40.39	300	150	45
400.40.50	400.40.51	400.40.53	400.40.55	400.40.54	300	150	45
400.40.60	400.40.61	400.40.63	400.40.65	400.40.64	300	150	45

size	300x135x71mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	3.0255_3.3315 (Al 99,5% + AlMg1)
Screw	1.4305 (X8CrNiS18-9)
Handle	1.4301 (X5CrNi18-10)
Others	SILIKON (Poly(organo)siloxane)