








## Product data sheet REF 49.50.10

11.12.2025

General data			
article number	49.50.10	Medicontainer 290x290x100mm silver	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826051945
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	1,960	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.60 (catalogue plastic surgery ) 451.02.61 (brochure Medicontainer System ) 453.03.01 (Blepharoplasty - Instruments for blepharoplasty Brochure)		
Complete picture		Detail images	
  			
size	290x290x100mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
Storage conditions			
No special requirements			
Materials			
	<b>Material number (DIN designation)</b>		
<b>Main part</b>	3.0255_3.3535_3.3315 (Al 99,5% + AlMg3 + AlMg1)		
<b>Handle</b>	1.4301 (X5CrNi18-10)		
<b>Others</b>	SILIKON (Poly(organo)siloxane)		