








Product data sheet REF 49.50.17

11.12.2025

General data			
article number	49.50.17	Medicontainer 290x290x170mm silver	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826051976
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	2,300	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.60 (catalogue plastic surgery) 451.02.61 (brochure Medicontainer System)		
Complete picture		Detail images	
  			
size	290x290x170mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.0255_3.3535_3.3315 (Al 99,5% + AlMg3 + AlMg1)		
Handle	1.4301 (X5CrNi18-10)		
Others	SILIKON (Poly(organo)siloxane)		