



## Product data sheet REF 49.50.22

11.12.2025

General data			
article number	49.50.22	Medicontainer 290x290x220mm silver	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826051983
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	2,520	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	



Medicon media  
450.00.60 (catalogue plastic surgery )  
451.02.61 (brochure Medicontainer System )

## Complete picture

## Detail images



Article No.	Article No.	Article No.	Article No.	Article No.	Article No.	Article No.	Article No.	Article No.	Article No.
49.50.11	49.50.12	49.50.13	49.50.14	49.50.15	49.50.16	49.50.17	49.50.18	49.50.19	49.50.20
49.50.21	49.50.22	49.50.23	49.50.24	49.50.25	49.50.26	49.50.27	49.50.28	49.50.29	49.50.30
49.50.31	49.50.32	49.50.33	49.50.34	49.50.35	49.50.36	49.50.37	49.50.38	49.50.39	49.50.40
49.50.41	49.50.42	49.50.43	49.50.44	49.50.45	49.50.46	49.50.47	49.50.48	49.50.49	49.50.50

size	290x290x220mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
------	---------------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	3.0255_3.3535_3.3315 (Al 99,5% + AlMg3 + AlMg1)
Handle	1.4301 (X5CrNi18-10)
Others	SILIKON (Poly(organo)siloxane)