



Product data sheet REF 49.50.27

11.12.2025


General data			
article number	49.50.27	Medicontainer 290x290x270mm laundry silver	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826051990
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	2,800	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	



Medicon media	451.02.61 (brochure Medicontainer System)
---------------	--

Complete picture

Detail images



size	290x290x270mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
------	---------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	3.0255_3.3535_3.3315 (Al 99,5% + AlMg3 + AlMg1)
Handle	1.4301 (X5CrNi18-10)
Others	SILIKON (Poly(organo)siloxane)