






Product data sheet REF 49.60.27

11.12.2025

| General data | | | |
|---|---|--|--|
| article number | 49.60.27 | Medicontainer 590x290x270mm laundry silver | |
| Unit | piece | Customs code (HS code) | 76169990 |
| GMDN | 13730 | GTIN | 04046826053093 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 13-730 |
| Classification | 1 (EU MDR) | Base UDI-DI | 4046826;N003;04;XB |
| Net weight (kg) | 4,780 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 451.02.61 (brochure Medicontainer System) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 590x290x270mm | Coating | Artificially created oxide layer (interference colors) on the surface of aluminum products |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.0255_3.3535_3.3315 (Al 99,5% + AlMg3 + AlMg1) | | |
| Handle | 1.4301 (X5CrNi18-10) | | |
| Others | SILIKON (Poly(organo)siloxane) | | |