








## Product data sheet REF 49.62.10

11.12.2025

General data			
article number	49.62.10	Medicontainer 590x290x110mm silver double lid	
Unit	piece	Customs code (HS code)	76169990
GMDN	13730	GTIN	04046826053550
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	4,800	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	451.02.61 (brochure Medicontainer System )		
Complete picture		Detail images	
  			
size	590x290x110mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.0255_3.3535_3.3315 (Al 99,5% + AlMg3 + AlMg1)		
Handle	1.4301 (X5CrNi18-10)		
Others	SILIKON (Poly(organo)siloxane)		