







Product data sheet REF 49.66.15

11.12.2025

General data																															
article number	49.66.15	bottom non perf 590x290x135mm																													
Unit	piece	Customs code (HS code)	76169990																												
GMDN	13730	GTIN	04046826053963																												
Risk category	1 (EU MDR)	UMDNS-Code	13-730																												
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB																												
Net weight (kg)	2,160	Operating instructions	GATDN00X																												
Latex free	yes	Preparation instructions																													
 																															
Medicon media	451.02.61 (brochure Medicontainer System)																														
Complete picture		Detail images																													
  <table border="1" data-bbox="1244 896 1452 1008"> <thead> <tr> <th>No.</th> <th>a (mm)</th> <th>b (mm)</th> <th>c (mm)</th> </tr> </thead> <tbody> <tr> <td>49.66.10</td> <td>590</td> <td>290</td> <td>85</td> </tr> <tr> <td>49.66.13</td> <td>590</td> <td>290</td> <td>120</td> </tr> <tr> <td>49.66.15</td> <td>590</td> <td>290</td> <td>135</td> </tr> <tr> <td>49.66.17</td> <td>590</td> <td>290</td> <td>155</td> </tr> <tr> <td>49.66.22</td> <td>590</td> <td>290</td> <td>205</td> </tr> <tr> <td>49.66.27</td> <td>590</td> <td>290</td> <td>255</td> </tr> </tbody> </table>				No.	a (mm)	b (mm)	c (mm)	49.66.10	590	290	85	49.66.13	590	290	120	49.66.15	590	290	135	49.66.17	590	290	155	49.66.22	590	290	205	49.66.27	590	290	255
No.	a (mm)	b (mm)	c (mm)																												
49.66.10	590	290	85																												
49.66.13	590	290	120																												
49.66.15	590	290	135																												
49.66.17	590	290	155																												
49.66.22	590	290	205																												
49.66.27	590	290	255																												
size	590x290x135mm	Coating	Artificially created oxide layer (interference colors) on the surface of aluminum products																												
Storage conditions																															
No special requirements																															
Materials																															
	Material number (DIN designation)																														
Main part	3.0255 (Al 99,5%)																														
Handle	1.4301 (X5CrNi18-10)																														