












Product data sheet REF 49.91.05

11.12.2025

| General data | | | |
|---|---|---|--|
| article number | 49.91.05 | paper filter w indicator 225x225mm neutral P/500 | |
| Unit | box/500 | Customs code (HS code) | 48054000 |
| GMDN | 13730 | GTIN | 04046826425425 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 13-734 |
| Classification | 16 (EU MDR) | Base UDI-DI | 4046826;N003;04;XB |
| Net weight (kg) | 1,678 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|     | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
| | | | |
| size | 225x225mm | Coating | |
| Storage conditions | | | |
| 10 °C |  30 °C | 10 % |  60 % |
| | |  |  |
|  | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| | | | |