












## Product data sheet REF 49.91.20

11.12.2025

General data			
article number	49.91.20	indicator label 63.5x36mm p/1000 adhesive	
Unit	box/1000	Customs code (HS code)	48025890
GMDN	13730	GTIN	04046826054250
Risk category	KEIN MP (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,550	Operating instructions	
Latex free	yes	Preparation instructions	
  			
Medicon media	451.02.61 (brochure Medicontainer System )		
Complete picture		Detail images	
			
size	63.5x36mm	Coating	
Storage conditions			
10 °C	30 °C	10 %	60 %
			
			
			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		