













Product data sheet REF 49.91.23

11.12.2025

General data			
article number	49.91.23	indicator label 50x18mm p/1000 adhesive	
Unit	box/1000	Customs code (HS code)	48025890
GMDN	13730	GTIN	04046826054281
Risk category	KEIN MP (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,406	Operating instructions	
Latex free	yes	Preparation instructions	
  			
Medicon media	451.02.61 (brochure Medicontainer System) 451.02.65 (brochure Medicontainer 2.0 Alu Container f Sterilization)		
Complete picture		Detail images	
			
size	50x18.0mm	Coating	
Storage conditions			
10 °C	 30 °C	 10 %	 60 %
			
			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		