














Product data sheet REF 49.91.31

11.12.2025

| General data | | | |
|---|--|---|---|
| article number | 49.91.31 | paper filter w indicator 265x120mm p/500 | |
| Unit | box/500 | Customs code (HS code) | 48054000 |
| GMDN | 13730 | GTIN | 04046826054304 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 13-734 |
| Classification | 16 (EU MDR) | Base UDI-DI | 4046826;N003;04;XB |
| Net weight (kg) | 2,200 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|     | | | |
| Medicon media | 451.02.61 (brochure Medicontainer System) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 265x120mm | Coating | |
| Storage conditions | | | |
| 10 °C | 30 °C | 10 % | 60 % |
|  | |  | |
| | |  |  |
| | |  | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | PAPIER KREPP WEISS 60G/M ² (Sterilisationspapier krepp weiß, 60g/m ²) | | |
| Others | WASSERDAMPF_INDIKATOR (Wasserdampf-Indikator) | | |