















Product data sheet REF 49.91.33

11.12.2025

General data			
article number	49.91.33	paper filter w indicator 240x120mm p/500	
Unit	box/500	Customs code (HS code)	48054000
GMDN	13730	GTIN	04046826054328
Risk category	1 (EU MDR)	UMDNS-Code	13-734
Classification	16 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	2,150	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
   			
Medicon media	451.02.61 (brochure Medicontainer System)		
Complete picture		Detail images	
			
size	240x120mm	Coating	
Storage conditions			
10 °C	30 °C	10 %	60 %
			
			
			
Materials			
	Material number (DIN designation)		
Main part	PAPIER KREPP WEISS 60G/M ² (Sterilisationspapier krepp weiß, 60g/m ²)		
Others	WASSERDAMPF_INDIKATOR (Wasserdampf-Indikator)		