



Product data sheet REF 49.91.67

11.12.2025

| General data | | | |
|-----------------|------------------|------------------------------|----------------|
| article number | 49.91.67 | ident label 50x17.5mm orange | |
| Unit | piece | Customs code (HS code) | 76169990 |
| GMDN | 60909 | GTIN | 04046826054489 |
| Risk category | KEIN MP (EU MDD) | UMDNS-Code | |
| Classification | NA (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,002 | Operating instructions | |
| Latex free | yes | Preparation instructions | MPU01 |



| | |
|---------------|--|
| Medicon media | 451.02.61 (brochure Medicontainer System) 451.02.65 (brochure Medicontainer 2.0 Alu Container f Sterilization) |
|---------------|--|

Complete picture

Detail images



| Farbe | ohne Gravur | mit Gravur |
|--------|-------------------|----------------|
| Color | without engraving | with engraving |
| Color | sin grabaci3n | con grabaci3n |
| Colour | senza gravure | sen gravure |
| Colore | senza incisione | con incisione |
| | 50 x 17,5 mm | |
| ■ | 49.91.61 | 49.92.61 |
| ■ | 49.91.62 | 49.92.62 |
| ■ | 49.91.63 | 49.92.63 |
| ■ | 49.91.64 | 49.92.64 |
| ■ | 49.91.65 | 49.92.65 |
| ■ | 49.91.66 | 49.92.66 |
| ■ | 49.91.67 | 49.92.67 |
| ■ | 49.91.68 | 49.92.68 |

| | | | |
|------|-----------|---------|--|
| size | 50x17.5mm | Coating | |
|------|-----------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 3.0255 (Al 99,5%) |