






Product data sheet REF 49.97.70

11.12.2025

| General data | | | |
|---|--|----------------------------|--------------------|
| article number | 49.97.70 | silicone mat 3mm 460x240mm | |
| Unit | piece | Customs code (HS code) | 39269097 |
| GMDN | 13730 | GTIN | 04046826055387 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 12-142 |
| Classification | 1 (EU MDR) | Base UDI-DI | 4046826;N003;02;X5 |
| Net weight (kg) | 0,551 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 451.02.61 (brochure Medicontainer System) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 3mm 460x240mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | SILIKON (Poly(organo)siloxane) | | |