



## Product data sheet REF 50.17.73

11.12.2025

General data			
article number	50.17.73	Doyen AT intest clamp 23cm cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826056209
Risk category	1r (EU MDR)	UMDNS-Code	10-871
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,055	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1



Medicon media	450.00.03 (general catalogue the art of surgery)
---------------	--

## Complete picture

## Detail images



size	23cm	Coating	
------	------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)