







## Product data sheet REF 50.23.73

11.12.2025

| General data   |  |                                 |                         |
|--|--|---------------------------------|-------------------------|
| article number   | 50.23.73   | Kocher AT intest clamp 23cm cvd |                         |
| Unit   | piece  | Customs code (HS code)          | 90189084                |
| GMDN   | 62468  | GTIN                            | 04046826056391          |
| Risk category  | 1r (EU MDR)                                      | UMDNS-Code                      | 10-871                  |
| Classification   | 6 (EU MDR)                                       | Base UDI-DI                     | 4046826;B;R002;04;ah;4K |
| Net weight (kg)  | 0,056  | Operating instructions          | GATDR00X                |
| Latex free   | yes  | Preparation instructions        | MPK23-1                 |
|      |  |                                 |                         |
| Medicon media  | 450.00.03 (general catalogue the art of surgery) |                                 |                         |
| Complete picture   |  | Detail images                   |                         |
|   |  |                                 |                         |
| size   | 23cm   | Coating                         |                         |
| Storage conditions   |  |                                 |                         |
| No special requirements  |  |                                 |                         |
| Materials  |  |                                 |                         |
|  | Material number (DIN designation)                |                                 |                         |
| Main part  | 1.4021 (X20Cr13)                                 |                                 |                         |