







## Product data sheet REF 50.23.78

11.12.2025

General data			
article number	50.23.78	Kocher AT intest clamp 28cm cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826056414
Risk category	1r (EU MDR)	UMDNS-Code	10-871
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,066	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
 			
size	28cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		