



Product data sheet REF 50.93.72

11.12.2025

| General data | | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| article number | 50.93.72 | plastic cone only f 50.93.65/70/78 | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 35157 | GTIN | 04046826057879 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 12-335 |
| Classification | 1 (EU MDR) | Base UDI-DI | 4046826;N001;06;WT |
| Net weight (kg) | 0,005 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| Storage conditions | | | |
|  140 °C | |   | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | POM (Polyoxymethylen_(Polyacetal/_Polyformaldehyd)) | | |