




## Product data sheet REF 50.93.75

11.12.2025

General data			
article number	50.93.75	lig ring latex-free p/100 f 50.93.65/68/70/78	
Unit	box/100	Customs code (HS code)	40169997
GMDN	35157	GTIN	04046826057909
Risk category	2a (EU MDD)	UMDNS-Code	12-335
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,013	Operating instructions	G509375
Latex free	yes	Preparation instructions	MPU02
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
Storage conditions			
   			
Materials			
	Material number (DIN designation)		
Main part	IR_617 (Polymer)		