



## Product data sheet REF 50.95.02

11.12.2025

| General data  |  |                              |                         |
|---|--|------------------------------|-------------------------|
| article number  | 50.95.02   | cotton a sponge carrier 28cm |                         |
| Unit  | piece  | Customs code (HS code)       | 90189084                |
| GMDN  | 32873  | GTIN                         | 04046826058050          |
| Risk category   | 1r (EU MDR)                                      | UMDNS-Code                   | 13-698                  |
| Classification  | 6 (EU MDR)                                       | Base UDI-DI                  | 4046826;B;R002;01;ah;3J |
| Net weight (kg)   | 0,030  | Operating instructions       | GATDR00X                |
| Latex free  | yes  | Preparation instructions     |                         |
|   |  |                              |                         |
| Medicon media   | 450.00.03 (general catalogue the art of surgery) |                              |                         |
| Complete picture  |  | Detail images                |                         |
|    |  |                              |                         |
| size  | 28cm   | Coating                      |                         |
| Storage conditions  |  |                              |                         |
| No special requirements   |  |                              |                         |
| Materials   |  |                              |                         |
|   | Material number (DIN designation)                |                              |                         |
| Main part   | 1.4310 (X10CrNi18-8)                             |                              |                         |
| Handle  | 1.4310 (X10CrNi18-8)                             |                              |                         |
| Others  | 1.4305 (X8CrNiS18-9)                             |                              |                         |