







## Product data sheet REF 51.13.05

11.12.2025

General data			
article number	51.13.05	endoscope Ø 4.0mm 300mm 12°	
Unit	piece	Customs code (HS code)	90189020
GMDN	36654	GTIN	04046826059620
Risk category	2a (EU MDD)	UMDNS-Code	14-289
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,200	Operating instructions	G776000
Latex free	yes	Preparation instructions	0
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
  			
size	4.0mm 300mm 12°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	PEEK (Polyetheretherketon )		
Others	SAPHIRGLAS (Glas aus Saphirkristallen)		