








Product data sheet REF 51.57.98

11.12.2025

General data			
article number	51.57.98	Gil Vernet retr 18mm 24cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826060510
Risk category	1r (EU MDR)	UMDNS-Code	12-028
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,051	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
  			
size	18mm 24cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		