






Product data sheet REF 51.75.30

11.12.2025

| General data | | | |
|---|--|---|--------------------|
| article number | 51.75.30 | Guyon catheter introducing instr cvd 37cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 38749 | GTIN | 04046826063016 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 10-762 |
| Classification | 5 (EU MDR) | Base UDI-DI | 4046826;N004;02;XG |
| Net weight (kg) | 0,025 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 37cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |