






## Product data sheet REF 52.25.51

11.12.2025

| General data  |  |                                   |                    |
|---|--|-----------------------------------|--------------------|
| article number  | 52.25.51   | Matzenauer vag retr 100x10mm 22cm |                    |
| Unit  | piece  | Customs code (HS code)            | 90189084           |
| GMDN  | 35352  | GTIN                              | 04046826064761     |
| Risk category   | 1 (EU MDR)                                       | UMDNS-Code                        | 13-666             |
| Classification  | 5 (EU MDR)                                       | Base UDI-DI                       | 4046826;N004;05;XR |
| Net weight (kg)   | 0,130  | Operating instructions            | GATDN00X           |
| Latex free  | yes  | Preparation instructions          |                    |
|   |  |                                   |                    |
| Medicon media   | 450.00.03 (general catalogue the art of surgery) |                                   |                    |
| Complete picture  |  | Detail images                     |                    |
|   |  |                                   |                    |
| size  | 100x10mm 22cm                                    | Coating                           |                    |
| Storage conditions  |  |                                   |                    |
| No special requirements   |  |                                   |                    |
| Materials   |  |                                   |                    |
|   | Material number (DIN designation)                |                                   |                    |
| Main part   | 1.4301 (X5CrNi18-10)                             |                                   |                    |