






Product data sheet REF 52.30.21

11.12.2025

| General data | | | |
|---|--|----------------------------------|--------------------|
| article number | 52.30.21 | Hegar uter dilator 21.0mm 18.5cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 11267 | GTIN | 04046826066017 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 11-266 |
| Classification | 5 (EU MDR) | Base UDI-DI | 4046826;N004;02;XG |
| Net weight (kg) | 0,090 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 21.0mm 18.5cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |