






## Product data sheet REF 52.31.03

11.12.2025

General data			
article number	52.31.03	Hegar uter dilator 3.5mm 18.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	11267	GTIN	04046826066147
Risk category	1 (EU MDR)	UMDNS-Code	11-266
Classification	5 (EU MDR)	Base UDI-DI	4046826;N004;02;XG
Net weight (kg)	0,015	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	3.5mm 18.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		