






Product data sheet REF 52.35.03

11.12.2025

General data			
article number	52.35.03	Hegar uter dilator d/e 3+4mm 20cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	11267	GTIN	04046826066765
Risk category	1 (EU MDR)	UMDNS-Code	11-266
Classification	5 (EU MDR)	Base UDI-DI	4046826;N004;02;XG
Net weight (kg)	0,017	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	3+4mm 20cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		