






Product data sheet REF 52.50.35

11.12.2025

| General data | | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------|------------------------------------------------|-------------------------|
| article number | 52.50.35 | Sommer uter elevating forceps 24cm cvd lateral | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 62468 | GTIN | 04046826067656 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 11-797 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R002;02;ah;3V |
| Net weight (kg) | 0,147 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 24cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |