






Product data sheet REF 52.71.01

11.12.2025

| General data | | | |
|---|--|---------------------------------------|-------------------------|
| article number | 52.71.01 | Schroeder uter scoop rigid 30cm fig 1 | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 46870 | GTIN | 04046826068509 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-508 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;02;ai;4K |
| Net weight (kg) | 0,076 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 30cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Handle | 1.4021 (X20Cr13) 1.4301 (X5CrNi18-10) | | |