



Product data sheet REF 52.73.08

11.12.2025

| General data | | | |
|-----------------|-------------|-------------------------------|-------------------------|
| article number | 52.73.08 | Gourdet uter curette 8mm 28cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 32594 | GTIN | 04046826068707 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 11-084 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;02;ai;4K |
| Net weight (kg) | 0,050 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK12-1 |



| | |
|---------------|--|
| Medicon media | 450.00.03 (general catalogue the art of surgery) |
|---------------|--|

Complete picture

Detail images



| | | | |
|------|----------|---------|--|
| size | 8mm 28cm | Coating | |
|------|----------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|---------------------|-----------------------------------|
| Main part Handle | 1.4021 (X20Cr13) |
| | 1.4301 (X5CrNi18-10) |