








## Product data sheet REF 52.80.03

11.12.2025

| General data  |  |                                       |                         |
|---|--|---------------------------------------|-------------------------|
| article number  | 52.80.03   | Sims uter curette rigid sh 26cm fig 3 |                         |
| Unit  | piece  | Customs code (HS code)                | 90189084                |
| GMDN  | 32594  | GTIN                                  | 04046826068868          |
| Risk category   | 1r (EU MDR)                                      | UMDNS-Code                            | 11-084                  |
| Classification  | 6 (EU MDR)                                       | Base UDI-DI                           | 4046826;B;R003;02;ai;4K |
| Net weight (kg)   | 0,083  | Operating instructions                | GATDR00X                |
| Latex free  | yes  | Preparation instructions              | MPK12-1                 |
|     |  |                                       |                         |
| Medicon media   | 450.00.03 (general catalogue the art of surgery) |                                       |                         |
| Complete picture  |  | Detail images                         |                         |
|    |  |                                       |                         |
| size  | 26cm   | Coating                               |                         |
| Storage conditions  |  |                                       |                         |
| No special requirements   |  |                                       |                         |
| Materials   |  |                                       |                         |
|   | Material number (DIN designation)                |                                       |                         |
| Main part   | 1.4021 (X20Cr13)                                 |                                       |                         |