








Product data sheet REF 52.80.05

11.12.2025

| General data | | | |
|---|--|---------------------------------------|-------------------------|
| article number | 52.80.05 | Sims uter curette rigid sh 26cm fig 5 | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 32594 | GTIN | 04046826068882 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 11-084 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;02;ai;4K |
| Net weight (kg) | 0,092 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK12-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 26cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |