








Product data sheet REF 53.70.10

11.12.2025

General data			
article number	53.70.10	Recamier placenta curette sh 40mm 35cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32594	GTIN	04046826070991
Risk category	1r (EU MDR)	UMDNS-Code	11-084
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;02;ai;4K
Net weight (kg)	0,088	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
  			
size	40mm 35cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Handle	1.4301 (X5CrNi18-10)		