





Product data sheet REF 53.80.04

11.12.2025

General data			
article number	53.80.04	diagnostic sheath Ø 5.5mm w 1 stopcock	
Unit	piece	Customs code (HS code)	90189084
GMDN	37086	GTIN	04046826071264
Risk category	2a (EU MDD)	UMDNS-Code	12-081
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,200	Operating instructions	
Latex free	yes	Preparation instructions	MPK03
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	Ø 5.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		