



Product data sheet REF 54.35.15

11.12.2025

General data			
article number	54.35.15	lung retr 132mm 26cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826072360
Risk category	1r (EU MDR)	UMDNS-Code	13-649
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,180	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.53 (catalogue cardio)		
Complete picture		Detail images	
			
size	132mm 26cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Handle	1.4301 (X5CrNi18-10)		